


Alpha Non Slip Cushion

Prevents
sliding and
gives
support



The Alpha Non Slip Cushion is placed underneath the thighs of a client seated in bed. A gentle and comfortable supporting and securing function is achieved through the non-slip properties of the surface.

It significantly prevents the client from sliding down in the bed while sitting.

Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Non Slip Cushion	PM-8050	40 x 50 x 8 cm 15,7" x 19,6" x 3,1"	Polyethylen, latex free	200kg / 440lb

    = 180 kg / 396lb


PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.

Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.