

# Alpha Sliding Mat


Non-static  
Low friction  
Soft



PM-2015



PM-2010




PM-1010

The Alpha Sliding Mats reduce the frictional resistance while moving the client. They are particularly helpful in easing client transfers in bed, into the shower trolley, the (wheel) chair and more.

Depending on the client's functional capacities, the Alpha Sliding Mat also supports the chest and pelvis for positioning and turning onto the side in a lying position.

## Technical specification



PM-1210



PM-1420


Product Name	Art. nr.	Dimensions	Material
Alpha Sliding Mat	PM-1010	190 x 60 cm 74,8" x 23,6"	PES/Carbon, latex free
Alpha Sliding Mat Medium	PM-1210	100 x 70 cm 39,4" x 27,5"	PES/Carbon, latex free
Alpha Sliding Mat Cotton Medium	PM-1220	100 x 70 cm 39,4" x 27,5"	Cotton/natural ticking, latex free
Alpha Sliding Mat Large	PM-1215	110 x 75 cm 43,3" x 29,5"	PES/Carbon, latex free
Alpha Sliding Mat Cotton	PM-1420	140 x 90 cm 55,1" x 35,4"	Cotton/natural ticking, latex free
Alpha Sliding Mat	PM-2010	50 x 60 cm 19,7" x 23,6"	PES/Carbon, latex free
Alpha Sliding Mat Cotton	PM-2020	50 x 60 cm 19,7" x 23,6"	Cotton/natural ticking, latex free
Alpha Sliding Mat 2 pieces	PM-2015	25 x 60 cm 9,8" x 23,6"	PES/Carbon, latex free

**PETERMANN**  
part of the Human Care Group

These products are distributed by Human Care AB for Petermann GmbH.



CE MDI = 200 kg / 440lb



Visit [www.humancaregroup.com](http://www.humancaregroup.com) for the latest product information and contact details to national offices and partners.