




OnlyU Disposable Sliding Foil

OnlyU Wall Mount & Rack



The OnlyU Disposable Sliding Foil is a tubular sheet which is skin-compatible, tear-resistant, latex free, biocompatible and food-safe.

This handling aid is used in clinics and hospitals for all the various patient position changes and transfers. The disposability of the OnlyU Disposable Sliding Foil means a rapid and particularly hygienic transfer for clients. The sliding foil can be cut into any desired length and be used just like a sliding mat.



Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
OnlyU Disposable Sliding Foil	PM-060.60-100	0,6 x 100 m 24" x 39,3"	Polyethylene (LDPE), latex free	200 kg 440 lb
OnlyU Disposable Sliding Foil incl. dispenser box & foil cutter	PM-060.65-100	0,6 x 100 m 24" x 39,3"	Polyethylene (LDPE), latex free	200 kg 440 lb
OnlyU Disposable Sliding Foil	PM-060.70-100	0,7 x 100 m 27,5" x 39,3"	Polyethylene (LDPE), latex free	200 kg 440 lb
Wall Rack for Sliding Foils	PMF-0200	46 x 70 x 18 cm / 18,1" x 27,5" x 7"	Stainless Steel	
OnlyU Wall mount for roll width 60 cm	PM-060.10			
OnlyU Wall mount for roll width 70 cm	PM-060.30			



PETERMANN
part of the Human Care Group

These products are distributed by Human Care AB for Petermann GmbH.



CE (MD) = 200 kg / 440 lb



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.