



This Alpha Turning Sheet with hidden handles was developed specially for clients who are particularly sensitive to pain. The processing of breathable textile products enables the client to remain on the Turning Sheet on an ongoing basis. The hidden handles and loops of the turn sheet make it possible for the carer to move the client freely and gently to the desired care situation in the bed. The Turning Sheet can be attached to the head of the bed or laterally to the slatted frame with the attachment straps provided.

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Turn Sheet with hidden handles	Cotton, natural thickening, latex free	Cream	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-2098	200 x 100	1650 g	(01)04250355320988(10)LOT

Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 90°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

