



The Alpha Non Slip Mat prevents the feet from slipping. When placed underneath the client's feet in bed, the non-slip mat helps to prevent slipping while repositioning or raising the patient. A non-slip mat increases the frictional resistance by 100%. It also makes it possible for the client moving toward the top of the bed by pressing with their feet.



Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Non Slip Mat	Polymeric Compound	Red	N/A		5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-7010	35 x 25	250 g	

Hygiene and maintenance

The product can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning. When necessary, we recommend wipe disinfection with the product Trionic® or comparable products. The use of bleaching agents or alcohol can permanently impair the properties of the product.



Made in GB