



BM-1600.3 + BM-1600.4

PM-7010



The Turn Sheet is a four-piece aid for reducing the friction forces during patient repositioning and movements in the bed. It consists of a Top Sheet a Fitted Sheet and two Non-Slip Mats. The Turn Sheet can remain in the bed permanently.

The Turn Sheet is very versatile. The sliding surface of the Fitted Sheet is finished with extremely low-friction silk. The fact that the silk layers of the Top Sheet and Fitted Sheet lie against one another

maximises the reduction of friction forces, producing the same effects as a sliding mat. Placing the Non-Slip Mats between the upper and lower sliding planes reduces slippage in this area.



Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Blue Motion Turn Sheet	Cotton/Polyamide latex free	anthracite/blue	180 kg	2 years	5 °C – 40 °C protect against sunlight
Blue Motion Non Slip Mat	Silicone rubber	Blue			5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
BM-1600 Set:			
BM-1600.3	200 x 140	675 g	(01)04250355316035(10)LOT
BM-1600.4	200 x 90	765 g	(01)04250355316042(10)LOT
PM-7010 2 pcs.	35 x 25	250 g	(01)04250355370105(10)LOT

Hygiene and maintenance

The Blue Motion Slide Sheet doesn't require any special maintenance. It can be machine washed at 60°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

BM-1600:



Made in Germany

The Non Slip Mat can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning. When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.

Non Slip Mat:



Made in UK